

## Rave ODM Adapter Training Workshop

### Instructor Led Training Course Outline

**Course Description:** Medidata’s ODM Adapter is a Rave Web Service capability which provides users with transactional ODM extracts of clinical transaction records for fields, field workflow status, and field queries for a given study. This solution enables partners and sponsors to load Rave EDC data into their enterprise warehousing system. This workshop will take users through a series of exercises that will teach them to effectively utilize Medidata’s ODM Adapter.

Approximate Duration: 2 ½ day sessions (8 hours)

Time	Topics
<b>Session 1</b>	
15 mins	<p><b>Welcome</b></p> <p>Introductions</p> <p>Course Overview</p>
15 mins	<p><b>Introduction</b></p> <p>Introduction to RWS and ODM Adapter</p> <p>Vendor Extensions</p> <p>How Rave Objects Map to ODM</p>
1 hour 45 mins <sup>1</sup>	<p><b>Clinical Data</b></p> <p>Subject Listing Dataset</p> <p>Introduction to Clinical Audits Records Dataset</p> <p>How to Use Pagination When Calling the Clinical Audits Record Dataset</p> <p>Understanding the Clinical Audits Sequence of Events</p> <p>Understanding Clinical Audits Records Audit Sub-Categories</p> <p>Inclusions and Exclusions</p>

---

<sup>1</sup> A 15 minute break will be included at the trainer’s discretion.

30 min	<b>Exercise 1</b> Log Lines
30 min	<b>Exercise 2</b> Queries
30 min	<b>Exercise 3</b> Statuses
15 min	<b>Q&amp;A</b>
<b>Session 2</b>	
2.5 hrs <sup>2</sup>	<b>Metadata</b> Study Metadata Dataset Understanding How Multiple CRF Versions Affect ODM Adapter Understanding How CRF Migrations Affect ODM Adapter Version Folders Dataset Signatures Dataset Local Lab Analyte Ranges Dataset
30 min	<b>Exercise 4</b> Field Restrictions
15 min	<b>Admin Data</b> Sites Listing Dataset Users Listing Dataset
15 min	<b>Best Practices</b>
15 min	<b>Q&amp;A</b>

---

<sup>2</sup> A 15 minute break will be included at the trainer's discretion.